Poster Presentation for JALAS 2019

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Construction of experimental environment for Microbiome animal experiments using the bioBUBBLE Enclosure and the Individual Ventilation Cage

Purpose

The animal experiments with germ-free or gnotobiotic mice used to be operated with a vinyl isolator, but it takes much time to prepare the vinyl isolator. Recently, the rearing and maintenance of germ-free mice using individual ventilation cages has been promoted.

In order to apply the breeding technology to microbiome animal experiments, the bioBUBBLE enclosure and IVC are combined to construct a breeding environment.

The method comprises examining the possibility of keeping and simultaneously breeding germ-free mice and not germ-free mice simultaneously breeding and evaluating the ability to control microorganisms in the bioBUBBLE enclosure and IVC.

Assessment

Experiment 1: Confirming whether germ-free mice can be kept and maintained in a positive pressure-controlled conventional breeding room. (3 cages)

Experiment 2: Confirming whether germ-free mice and ASF mice can be reared in the same IVC rack. (2 cages for germ-free mice, 4 cages for ASF mice)

Experiment 3: In the BSL 2 lab room, where the microbiome experiment is conducted under negative pressure, for breeding germ-free mice by using the bioBUBBLE enclosure and safety cabinet.

Materials and Methods

Construction of experimental environment using the bioBUBBLE enclosure and IVC

1) The conventional breeding room under the positive pressure (Figure 1)

Install the bioBUBBLE enclosure with a working station and pass-through box. The IVC is made by Japan CLEA. (Figure 2). The operator must wear dust-free clothes in order to enter the bB enclosure, and must also change the cages in the working station.

2) BSL 2 lab room under negative pressure (Figure 3)

Assuming the environment of the microbiome animal experiment, install the bioBUBBLE enclosure in order to cover the safety cabinet. The IVC rack is made by Japan CLEA (Figure 2). The operator must wear dust-free clothes in order to enter the bB enclosure,

and must also change the cages in the safety cabinet.

Equipment sterilization method

One set of breeding equipment was prepared (Figure 4) for each IVC. The cage was placed in a double sterile bag with the lid of the cage shifted in order to improve sterilization conditions (Figure 5). The sterilization of breeding equipment was at 127 degree C for 40 minutes (Figure 6). In the sterilization of water, the heat-resistant container was filled with warm water, and the container was placed in a double sterilization bag with the lid of the container loosened (Figure 7). From the same preliminary study, sterilization conditions were set at 127 degree C for 90 minutes. The sterilized equipment was removed and stored in the pass-box, sprayed with MB-10 (500 ppm), and left overnight.

The method of introducing IVC to the working station

MB-10 was sprayed into the working station before work. The IVC in use was introduced into the working station by spraying MB-10 on the exterior (Figure 8). Due to the space of the working station, one loading was upto 3 cages. We sprayed MB-10 in the working station and let it stand for 5 minutes, and then carried in the sterilized IVC from the pass-box.

Cage replacement

The mice were moved using tweezers. The tweezers were placed in a 70% ethanol-containing tweezer stand and exchanged for each cage. The hands were disinfected using a hand spray containing 70% ethanol each time (Figure 9). The cage was changed once a week. In Experiment 2, the cage replacement was carried out for 3 cages of ASF mice first, next, cage replacement was performed in the order of germ-free, ASF, and germ-free (Figure 10). Tweezers and a hand spray were used for both germ-free and ASF mice, and the gloves were changed when operating the IVCs for germ-free mouse breeding.

Cleaning and Disinfection

After work, we cleaned and disinfected the working station, work bench, and floor with 70% ethanol (Figure 11). The outside of the bioBUBBLE enclosure was cleaned and disinfected using chlorinated water (Figure 12).

Germ-free inspection and ASF mice

The sterility test conforms to the test method of the Japan Society of Experimental animals once a week, and it was carried out in all the cages. For the ASF mice, we used the Altered Schaedler Flora (ASF) established gnotobiotic mice. ASF inspection was also performed as frequently as sterility testing.

Results, Consideration, Supervision

Experiment 1

Result: 38 weeks of breeding germ-free mice was confirmed, and maintenance was terminated.

Consideration: In microbiome animal experiments, the period from bacterial administration to evaluation requires approximately 1 month, and we believe that further long-term breeding is also possible.

Experiment 2

Result: Mice were bred for 29 weeks with germ-free IVC cages and ASF IVC cages mixed in the same rack, and maintenance was terminated.

Consideration: The results of maintaining germ-free mice indicate that they can be managed without cross-contamination between IVC cages. If it is conventional, we need to prepare the vinyl isolator for each bacteria to be evaluated, by using IVC cages; it was suggested that a significant reduction of working space and time could be achieved.

Experiment 3

Result: Maintenance of 14 weeks of germ-free mice breeding was confirmed and terminated.

Consideration: Based on the results of experiments 1 and 2, we proceeded to evaluate mice in the BSL 2 environment, where microbiome animal experiments will be conducted. In the BSL 2 environment, there is a need for negative pressure, and the mice in the safety cabinet are kept under negative pressure (Figure 13). So, we set up the bioBUBBLE enclosure to cover the safety cabinet; it was possible to maintain the germ-free mice.

Supervision

By using both the bioBUBBLE enclosure and the IVC cages, a new basic technology applicable to microbiome animal experiments has been developed. As a further step, we will proceed with the verification of reliability in actual microbiome animal experiments. At the same time, we aim to create the SOP and spread the technology, aiming at generalization of this technology.