


PU1515 HEPA Air Handler  
 Recirculates Within Larger  
 Negative Pressure Enclosure

Hard Room Ceiling  
 Must Have Min 3"  
 Clearance Above Exhaust

Magnetic Patch Openings  
 Along Front of Enclosure  
 (Flap Downward)

Negative/Front Patch:  
 26"W x 14"H  
 Positive Patch  
 22"W x 12"H

Entire Vinyl Side (All Sides)  
 Opens Making Internal  
 Equipment Accessible



Enclosure for Product and Personnel Protection
Enclosure OD: 58"W x 30"D x 52"H

Drawing Is Accepted By:
Signature:
Date:

<b>BIO BUBBLE</b>	
CONTROLLED ENVIRONMENTS   CUSTOM SOLUTIONS	
Phone: 970.224.4262 Fax: 970.224.2419 Email: sales@biobubble.com	
Version #:	02
Date:	14 Nov 2018
Scale:	No Scale
Drawn By:	Jack