RACK COVERS

bio**BUBBLE** Rack Covers are custom manufactured to fit any size rack, trolley, cart, or other vivarium fixture. Rack Covers are constructed of durable vinyl coated polyester fabric with Velcro®-type or zippered connections creating the longest lasting covers on the market. Use Rack Covers for the transport of both clean and dirty materials to maintain the highest possible cleanliness levels in the vivarium. Rack Covers are available in 8 different colors. Color-coding is useful for identification of components within the facility. Rack Covers easily slip onto your existing racks and can be secured around the lowest shelf (optional). Rack Covers are washable on an alkaline cycle at temperatures up to 180° F.

FEATURES & BENEFITS

- Custom manufactured to fit your racks, carts or trolleys.
- Heavy-duty fabrics create the most durable covers on the market.
- Washable in the cage / rack washer up to 180° F
- Rack Covers can be docked directly to bioBUBBLE
 Clean Rooms, isolators, and autoclave containers for completely isolated transfer of clean materials and animals. (optional)





HOW TO ORDER

You **must** measure your rack very accurately. The cover will slip down over your rack. Please follow these measuring procedures when ordering:

- **Step 1:** Lower the bottom shelf on your rack to its lowest position.
- Step 2: Measure the height from the bottom of the lowest shelf to the top of the top shelf.
- **Step 3:** Measure the width and depth, including *any protrusions*. Your Rack Cover will have approximately .5 1 IN (1.27 2.54 CM) of clearance on all the dimensions.
- Step 4: Write the dimensions on the drawing below.
- Step 5: Indicate the quantity for each color.

Standard	white	yellow	red	green
	black	royal blue		

Step 6: Choose Closure Type:

hook and loop (velcro-type)
zippered closures

- Step 7: Do you want hook and loop around the inside bottom of the cover to attach to the bottom shelf of your rack? Yes No (Circle one)
- **Step 8:** Return this form with your order.



